

Difco Middlebrook 7H10 Agar

1. Suspend 19 g of the powder in 900 mL of purified water containing 5 mL of glycerol. Mix thoroughly.
2. Autoclave at 121°C for 20 minutes.
3. Aseptically add 100 mL of Middlebrook OADC Enrichment to the medium when cooled to 50-55°C.
4. Pipette 35 ml per plate (petri dishes).

Difco Middlebrook 7H11 Agar

1. Suspend 21 g of the powder in 900 mL of purified water containing 5 mL of glycerol. Mix thoroughly.
2. Autoclave at 121°C for 20 minutes.
3. Aseptically add 100 mL of Middlebrook OADC Enrichment to the medium when cooled to 50-55°C.
4. Pipette 35 ml per plate (petri dishes).

Difco Middlebrook 7H9 Broth

1. Suspend 4,7 g of the powder in 900 mL of purified water containing 5 mL of Tween 80.
Mix thoroughly.
2. Autoclave at 121°C for 20 minutes.
3. Aseptically add 100 mL of Middlebrook ADC Enrichment to the medium when cooled to 45-50°C.