

# Stitching in ZEN Blue (Widefield)

1. Select **Processing** tab
2. Select **Stitching** under **Method** ( if it is not in your recently used list then it can be found under **Geometric** heading)
3. Select image to stitch (click on white arrow and select from the list)
4. Set output type in the **Parameters** tab to:
  - A. New Output
  - B. Fuse Tiles(This only needs to be done at the start of each session)
- 5 Press **Apply**

This will create a new image in the image list that is a fused version of the raw (titled IP.....)

Then just “save as”

NB it is worth keeping the original raw version as well as the fused stitched in case you want to re-use your settings (you can't re-use a fused file) and also re-stitch

