

NMR Service request form

Samples must be submitted in solution ready to run in externally clean NMR tubes of Precision grade (507 minimum) and appropriate length.

<p>Name:</p> <p>Supervisor:</p> <p>NOMAD login:</p> <p>e-mail:</p>	<p>Sample Reference:</p> <p>(MUST be clearly labelled on tube)</p> <p>Solvent:</p> <p>Experiments Required:</p> <p>Expected Outcome:</p>
<p>The sample must be a clear solution with no suspended particles etc.</p> <p>Structure:</p>	<p>Molecular Weight:</p> <p>Quantity:</p> <p>Stability:</p> <p>Toxicity and other hazards (remember we may need to clear up spillages):</p>